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Summary of the pre-analytical phase regarding urinary specimen collection, transportation and handling.

Description	Reference
Urinary specimen storage time	(Leber, 2016)
- Room temperature: max 2h - Refrigerator (2-8 °C): max 24h	
Replication limit	(Jorgensen et al., 2015; Sharp, 2009); BILULU Expert opinion; (Leber, 2016)
- One sample per sample type every 24h	
Rejection or acceptance of samples for diagnosis of UTI?	
- Acceptable sample types	
○ Clean catch-voided midstream urine	
○ Supra-pubic puncture	
○ Indwelling catheter	
○ Single sondage	
○ Urinary stoma*	
○ Pedibag	
- Suboptimal sample types	
○ Samples preserved > 2h at room temperature or > 24h at 2-8 °C	
○ Samples in leaking containers	
- Unacceptable sample types	
○ 24-h urine collections	
○ Catheter tip/bag	
○ First-voided urine fraction	
○ Samples contaminated with faeces	
Definition of pyuria	(Roggeman and Zaman, 2001; Manoni et al., 2013)
- Flow cytometry techniques	
○ Adults: ≥ 20-25 WBC/μL	
○ Pediatrics (0-2 years of age): ≥ 10-25 WBC/μL	
- Microscopy techniques:	
○ All patients: 10 WBC/μL	

* High contamination rate; Abbreviations: WBC: white blood cells.

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Summary of urinary specimen processing.

Description	
Volume of the inoculation?	- Routine samples: at least 1 μL - Detection of yeasts and invasive samples: 10 μL
Plates that have to be inoculated?	- Routine samples: combination of a blood-containing agar plate and a MacConkey or selective/ chromagar plate - Detection of yeasts: additional inoculation of a sabouraud or chromagar plate
How to incubate?	- Incubation temperature and atmosphere ○ 35-37 °C ○ Blood-containing agar: 5-10% CO ₂ -atmosphere ○ Selective plate/ chromagar: ambient temperature (35-37 °C) ○ Biplate: ambient temperature (35-37 °C), re-incubation in 5-10% CO ₂ -atmosphere - Duration of incubation: ○ Midstream of outpatients < 65 years: at least 18 h ○ Midstream of outpatients ≥ 65 years: at least 36 h ○ Inpatients: at least 36 h ○ Indwelling catheter: at least 36 h ○ Supra-pubic puncture: at least 42 h ○ Culture of yeasts: at least 42 h

*Invasive samples: supra-pubic puncture and single sondage specimens retained as invasive samples.

*Additional inoculation of a 1 ul. plate can facilitate the count.

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Classification of the different micro-organisms.

Category 1: urogenital skin flora	Viridans streptococci Commensal <i>Neisseriaceae</i> <i>Lactobacillus</i> species Coagulase-negative staphylococci other than <i>S. saprophyticus</i> <i>Corynebacterium</i> species other than <i>C. urealyticum</i> <i>Aerococcus</i> other than <i>A. urinae</i> and <i>A. sanguinicola</i>
Category 2: (common) uropathogens	Gram-negative bacilli <i>S. aureus</i> <i>S. saprophyticus</i> <i>S. pneumoniae</i> Large colony Beta-hemolytic streptococci* (excl. <i>Streptococcus anginosus</i>) Enterococci Yeasts (<i>C. albicans</i> , <i>C. glabrata</i>) <i>Actinotignum schaalii</i> <i>Aerococcus urinae</i> <i>Aerococcus sanguinicola</i> <i>Corynebacterium urealyticum</i> <i>Gardnerella vaginalis</i>
Category 3^b: rare or unusual uropathogens	

* *Streptococcus agalactiae* in women of childbearing age (15–50 years of age) and neonates should always be mentioned as isolate.

^b This group of micro-organisms is considered as uropathogen if "pure" culture or if 10 times higher concentrations of these pathogens have been grown.

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Scheme for processing and workup of urine cultures based on method of collection, and number of pathogens.

Type of collection	Culture		Extent of workup	
	Colony count	Extent of workup	Colony count	Extent of workup
Midstream Straight catheter Pedibag	$\geq 10^4$ CFU/mL:	ID and AB	1 isolate $\geq 10^5$ CFU/mL and 2 isolates $< 10^4$ CFU/mL	ID and AB of the predominant uropathogen ID of the other uropathogens No workup, ask new sample
	$< 10^4$ CFU/mL:	ID and ask control sample (AB if pyuria) ID and AB	Other cases:	
Indwelling catheter Supra-pubic sondage Nephrostomy	$\geq 10^4$ CFU/mL:	ID and AB	1 isolate $\geq 10^5$ CFU/mL and 2 isolates $< 10^4$ CFU/mL	ID and AB of the predominant uropathogen ID of the other uropathogens No workup, ask new sample
	$< 10^4$ CFU/mL:	ID and ask control sample	Other cases:	
Supra-pubic aspiration	Always ID and AB, except for yeasts (only ID). Contact clinician for clinical information*			

*In case of exceptional results, e.g. presence of urethral flora or ≥ 3 different types of micro-organisms. Abbreviations: ID: identification; AB: antimicrobial susceptibility testing; CFU: colony forming unit.

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Additional comments that can be added to a clinician report.

Additional comments in case of:	
Patients with an indwelling catheter:	Polymicrobial infections are possible in patients with an indwelling catheter. Contact the laboratory if further workup is desirable. Catheter-associated cystitis is treated with antibiotics if signs and symptoms of infection are present. Antibiotic therapy should be started after elimination of the catheter. If the catheter is needed and symptoms of infection are present, a change of catheter should be performed under appropriate antibiotic therapy.
A negative culture of a pedibag urinary sample: Isolation of <i>Actinotignum schaalii</i>	A negative culture has a highly negative predictive value for absence of infection Urogenital commensal, although probable cause of UTI in predisposed patients, especially patients aged older than 60 years or with underlying pathology. Antibiotic of choice is amoxicillin. The pathogen is intrinsically resistant to ciprofloxacin and co-trimoxazol (Kristiansen R et al., 2014)
Isolation of <i>Corynebacterium urealyticum</i>	This pathogen is associated with cystitis induced by stone formation and encrusted cystitis or pyelitis (e.g. in cases of kidney transplantation (Lopez-Medrano et al., 2008)). Treatment consists of antibiotic combination therapy, acidification of urine for several weeks and surgical removal of bladder stone if present.
Isolation of <i>Corynebacterium</i> spp. (other than <i>C. urealyticum</i>)	This group of bacteria belong to the urogenital skin flora, but could be a cause of UTI. Results need to be interpreted in a clinical context. The following <i>Corynebacterium</i> spp. are associated with UTI (Jorgensen et al., 2015): <ul style="list-style-type: none"> • <i>C. amycolatum</i> • <i>C. aurimucosum</i> • <i>C. glucuronolyticum</i> • <i>C. minutissimum</i> • <i>C. riegeii</i> • <i>Arthrobacter</i> spp. • <i>C. pseudogenitalium</i> • <i>C. striatum</i>
<i>Staphylococcus aureus</i>	<i>S. aureus</i> is not a classical uropathogen. Isolation could be due to post-operative infection? Post- uromanipulation? Hematogenous spreading? Presence of foreign body material? Further investigation is necessary (BILULU Expert opinion)
<i>Actinotignum</i> spp.	This group of bacteria belong to the urogenital skin flora, but are associated with complicated UTI (Jorgensen et al., 2015).
<i>Aerococcus urinae</i>	Mostly contaminant but possible uropathogen in predisposed patients (prostate hypertrophy, elderly patients). Nitrofurans and amoxicillin are therapeutic options in case of infection. Sensitivity to co-trimoxazole is contested. There are no interpretative criteria available for antimicrobial susceptibility testing (Jorgensen et al., 2015).
<i>Aerococcus sanguinicola</i>	Mostly contaminant, particularly associated with urosepsis. Amoxicillin is the first treatment option. Mostly, this pathogen is susceptible to cefuroxim, cefotaxim, erythromycin, tetracyclines and linezolid. There are no interpretative criteria available for antimicrobial susceptibility interpretation (Jorgensen et al., 2015).

Reference:

1: Oyaert M, Van Meensel B, Cartuyvels R, Frans J, Laffut W, Vandecandelaere P, DeBeenhouwer H; BILULU Study Group. Laboratory diagnosis of urinary tract infections: Towards a BILULU consensus guideline. J Microbiol Methods

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Overview of the different comments that can be presented on the laboratory report.

Observation culture	Report to the clinician
Only urogenital	Report: 'Urogenital flora' (+ number of CFU/mL)
If urogenital flora > uropathogens	Report: "Predominantly urogenital flora with possible presence of urinary pathogens. If clinically relevant, please send a control sample to the laboratory." (+ number CFU/mL)
Uropathogens ≥ urogenital flora	Report 'urogenital flora' (+ number of CFU/mL) and "uropathogens" (+ number of CFU/mL)
Only uropathogens (≤ 3 types)	Report all pathogens (+ number of CFU/mL)
Mixed flora with > 3 types of uropathogens	Report: "3 or more different types of pathogens, probably contamination. If clinically relevant, please send a control sample to the laboratory" (+ number of CFU/mL)